

New
Innovation

Drummond Portable Pipet-Aid® XL Pipet Controllers

The Ergonomic Design is Proven to Reduce Neck and Upper Arm Strain When Pipetting Under a Hood

- Longer, lightweight handle lowers your arm position to a more comfortable position
- New quiet, more efficient pump extends operation period between charges
- Can be charged while in use to eliminate downtime
- Three-speed settings provide more precise control for use in the pipet filler mode and the pipet dispenser mode

Pipetting under a hood doesn't need to be a pain in the neck. The new Drummond Portable Pipet-Aid XL has a longer handle to reduce the amount of arm lift required to perform the same operation as with a conventional pipettor. The ergonomic design includes an adjustable sliding hand rest and a removable stand which enables the unit to be set down without contaminating the pipet. The Portable Pipet-Aid XL also features a new quieter, more powerful pump, and a new EZ Grip nosepiece. The unit is supplied complete with power supply/charger, 4 extra filters, and a holster-type wall bracket.

DRUMMOND PORTABLE PIPET-AID® XL PIPET CONTROLLERS

Cat. No.	Description
4-000-105	Portable Pipet-Aid XL 110V
4-000-205	Portable Pipet-Aid XL 230V, European Plug
4-000-205-A	Portable Pipet-Aid XL 240V, Australia Plug
4-000-205-E	Portable Pipet-Aid XL 240V, UK Plug

Replacement Parts

Cat. No.	Description
4-000-002	Rubber Insert (2/pk)
4-000-051	TC Filters (5/pk)
4-000-057	TC Gaskets (pair)
4-000-036	Battery Pack (XL only)
4-000-065	XL Power Supply/Charger, 110V
4-000-079	EZ Grip Nosepiece



The Developers of the Original Pipet-Aid

For more information, call 800/523-7480
or visit our web site at www.drummondsci.com



DRUMMOND
SCIENTIFIC COMPANY

500 Parkway, Box 700
Broomall, PA 19008
Call Toll Free: 1-800/523-7480
1-610/353-0200
FAX: 1-610/353-6204
Web site: www.drummondsci.com

