



Drummond Scientific Company

500 PARKWAY, BROOMALL, PENNSYLVANIA 19008-4203

PHONE : 610-353-0200 FAX 610-353-6204

www.drummondsci.com

February 23, 2024

Certificate of Conformance & Sterility

Products listed conform to Drummond Scientific Company specifications and drawings, and have passed inspection and testing criteria for sterility. This item has been manufactured in accordance with FDA GMP's as Medical Device Materials of construction have passed USP Class VI-121 C. plastics test for biosafety

<u>Item</u>	<u>Part #</u>	<u>Lot Number</u>	<u>Exp</u>
Sterile Self-Locking Filters (5pk)	4-000-051-S	82064	2026-09-20
Method Of Sterilization	Gamma Radiation		
Hydrophilic Pore Size	.8 um PES		
Hydrophobic Pore Size	.8 um PTFE		
Filter Diameter	24mm		
Housing Material	Polypropylene		
Bubble Point	>1 barnq		

Alexis Piccirilli
Quality Assurance